MSE攪拌子 PAT. **NEW**

MSE Rotator (Patented)

(MSE=Multi-Stacked Elements)

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Faster mixing than ever!!
The best choice for mixing in R&D!

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The MSE rotator achieves uniform mixing in a short time with no generation of retention regions in a beaker.

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
The support and the stacked mixing elements are detachable to change the number of mixing elements depending on the size of a beaker.

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
A stirring bar, a typical tool to mix fluids in a small vessel such as a beaker, generates insufficient mixing regions of the fluid in the central portion and lengthens the mixing time. The MSE rotator solves this problem with its unique structure.

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

Mixing elements A (open at inner circumference) and mixing elements B (open at outer circumference) are stacked alternately. The mixing elements A and B are designed with overlapping holes allowing the fluid sucked into the hollow portion of the MSE rotator to flow radially in the stacked mixing elements and to discharge from the side openings.

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The following compares the mixing times in the decolorization test using a stirring bar with a length of 60 mm and MSE rotator with an outer diameter of 54 mm and 37 mm (Rotation speed: 200 rpm/ beaker capacity: 3L).

Remarks!: Mixing time will be varied by such as rotation speed of rotators or stirrers, beaker capacity, physical properties of liquid.

	開始時 [Start]	3秒後 [3s]	5秒後 [5s]	10秒後 [10s]	15秒後 [15s]	20秒後 [20s]	30秒後 [30s]	45秒後 [45s]	60秒後 [60s]
<p>MSE攪拌子(外径54mm) MSE rotator (O.D.: 54mm)</p>									
<p>MSE攪拌子(外径37mm) MSE rotator (O.D.: 37mm)</p>									
<p>棒状攪拌子(長さ60mm) Stirring bar (60mm length)</p>									

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●Our products are manufactured under appropriate quality control. However, ISEL will not be liable in any way for any damage incurred by forces of nature or by the customer that are caused by misuse of using the products. The contents of this catalog are subject to change without prior notice.